

# ELECTRONIC SOMATIC CELL COUNT

## Fossomatic™ 90 (Raw Commingled Cow, Sheep, Goat, Water Buffalo and Camel Milk) IMS #16

(Unless otherwise stated all tolerances  $\pm 5\%$ )

1. **Laboratory Requirements (see Cultural Procedures (CP), items 33 & 34)** \_\_\_\_\_
  - a. Un-preserved samples may be tested up to 72 hours after initial collection \_\_\_\_\_
  - b. Samples may be tested up to 7 days after initial collection if preserved with 0.02% 2-bromo-2-nitropropane-1,3-diol (Bronopol™) or 0.05% potassium dichromate ( $K_2Cr_2O_7$ ) \_\_\_\_\_
  
2. **Comparative Test with DMSCC**  
**[NOT required as a co-requisite for certification of analysts in laboratories purchasing standards from a CERTIFIED provider (item 17.b)]** \_\_\_\_\_
  - a. Analysts certified for DMSCC \_\_\_\_\_
  - b. Each analyst seeking certification for the ESCC test shall perform the comparative test \_\_\_\_\_
    1. Test 4 samples (100K-200K, 300K-500K, 600K-800K and 900K-1.2M) in triplicate for both DMSCC (three separate smears each) and ESCC \_\_\_\_\_
    2. Results must be evaluated by the FDA/LPET LEO or LEO and shown to be acceptable prior to official use of test in laboratory \_\_\_\_\_
    3. Copy of comparison and results in QC record (or easily accessible on file in the laboratory); kept for as long as analyst is certified \_\_\_\_\_
  - c. Required for laboratories preparing in house standards or using commercially prepared standards (items 17.a and c) and for those testing goat or camel milk \_\_\_\_\_

### APPARATUS

3. **See CP items 1-4** \_\_\_\_\_
4. **Fossomatic 90 Electronic Somatic Cell Counter** \_\_\_\_\_
5. **Water Bath** \_\_\_\_\_
  - a. Circulating and thermostatically controlled to 37-42°C \_\_\_\_\_

**6. Pipettor, fixed volume or electronic (\_\_\_\_\_)**

- a. Calibrate to deliver 500 µL milk (see CP item 6.e)
- b. Maintain records

**7. Pipettor Tips**

- a. Disposable, replace for each sample

**REAGENTS**

**8. Stock Dye Solution, 0.1% Ethidium Bromide (caution TOXIC, use gloves when handling and do not breathe dust)**

- a. Dissolve 1.0 g ethidium bromide (C<sub>21</sub>H<sub>20</sub>BrN<sub>3</sub>) in 1 L DI or MS water by heating to 40-60°C
- b. Store in light-proof, air-tight bottle no more than 60 days
- c. Lab Prep. Date: \_\_\_\_\_ Exp. Date: \_\_\_\_\_

**9. Stock Rinsing Solution, 1% Triton X-100**

- a. Dissolve 10 mL Triton X-100 in 1 L DI or MS water by heating to 60°C
- b. Store in air-tight container no more than 25 days
- c. Lab Prep. Date: \_\_\_\_\_ Exp. Date: \_\_\_\_\_

**10. Stock Buffer Solution, 0.025 M Potassium Hydrogen Phthalate**

- a. Dissolve 51.0 g KH phthalate and 13.75 g KOH in 10 L DI or MS water by heating to 40-60°C
- b. Add 150 mL 1% Triton X-100 (item 9), store less than 7 days in airtight container
- c. Lab Prep. Date: \_\_\_\_\_ Exp. Date: \_\_\_\_\_

**11. Ammonium Hydroxide (NH<sub>4</sub>OH) Solution, Reagent Grade, 25%**

**WORKING SOLUTIONS**

**12. Working Dye Solution/Zero Control (used within 7 days)**

- a. Dilute 26 mL stock dye solution (item 8.a) to 2.5 L with stock buffer solution (item 10.b)
- b. Lab Prep. Date: \_\_\_\_\_ Exp. Date: \_\_\_\_\_

**13. Working Rinsing Solution (used within 7 days)**

a. Add 10 mL stock rinsing solution (item 9) to 25 mL of 25% NH<sub>4</sub>OH and dilute to 10 L with DI or MS water

b. Lab Prep. Date: \_\_\_\_\_ Exp. Date: \_\_\_\_\_

**14. All solutions labeled with date prepared and expiration date**

**15. Optionally use manufacturer's reagent kits and instructions specific for each Instrument**

**START UP**

**16. Cell Counter**

a. Assure adequate volume of working solutions, not used beyond expiration date(s)

b. Turn on power and cycle at least six times

c. Blind count  $\leq 5$

d. Vacuum pressure setting minimum of -40 kPa

e. Dispenser filling time 4-5 sec

f. Intake filling time 3-4 sec

g. **IF ANY ABOVE PARAMETERS ARE OUT OF VARIANCE, CORRECT BEFORE PROCEEDING**

h. Maintain records on all parameters

**17. Milk Standards**

a. Commercially prepared: \_\_\_\_\_

Lot #: \_\_\_\_\_ Date Rcd: \_\_\_\_\_

1. Four standards in ranges 100K-200K, 300K-500K, 600K-800K and 900K-1.2M

2. Perform DMSCC in triplicate on each standard in set and average counts; maintain records

3. Perform DMSCC check in rotation by all certified analysts

4. Use standards within one week

Lab Exp. Date: \_\_\_\_\_

b. Certified provider: \_\_\_\_\_

Lot #: \_\_\_\_\_ Exp. Date: \_\_\_\_\_ Date Rcd: \_\_\_\_\_

1. Four standards in ranges 100K-200K, 300K-500K, 600K-800K and 900K-1.2M \_\_\_\_\_
2. Maintain copies of all provided DMSCC values \_\_\_\_\_
3. Measure and maintain records of temperature (0.0-7.5°C) of standards as received \_\_\_\_\_
4. Maintain copies of all correspondence regarding problems \_\_\_\_\_
5. Use standards by manufacturer's expiration date \_\_\_\_\_
6. Failed standards shall be verified with DMSCC \_\_\_\_\_
  - a. If no analysts certified for DMSCC then a new set of standards is required \_\_\_\_\_
  - b. Do not continue with official testing until the new standard(s) test(s) in range \_\_\_\_\_

c. Laboratory prepared (weekly) \_\_\_\_\_

1. Prepare from raw milk > 18 hours old preserved with 0.05% potassium dichromate ( $K_2Cr_2O_7$ ) \_\_\_\_\_
2. Or, preserve with 0.02% 2-bromo-2-nitropropane-1,3-diol (Bronopol™) \_\_\_\_\_
3. Standards **cannot** be preserved with formalin \_\_\_\_\_
4. Prepare 4 standards in ranges 100K-200K, 300K-500K, 600K-800K and 900K-1.2M; use within one week \_\_\_\_\_

Lab Prep. Date: \_\_\_\_\_ Exp. Date: \_\_\_\_\_

5. Perform DMSCC in triplicate on each standard prepared and average counts; maintain records \_\_\_\_\_
6. Perform DMSCC check in rotation by all certified analysts \_\_\_\_\_

d. Hourly Control Sample (instrument drift check) \_\_\_\_\_

1. Use one of the standards (items 17.a, b or c) in the 600-800K range, run in triplicate and determine average \_\_\_\_\_
2. Optionally, prepare sufficient control/sample 600-800K range, run in triplicate and determine average \_\_\_\_\_

## PROCEDURE

### 18. Testing Standards (each time instrument used)

- a. Heat standards to 37-42°C (using a temperature control) and test within 30 min of reaching temperature, use once and then discard; i.e., do not re-use
- b. Invert 10 times, pipet 500 µL into intake chamber within 3 min
- c. Run standards in triplicate and average the counts for each level; maintain records
- d. Each standard's average must be within 10% of the DMSCC (item 17) for that level, except within 15% for 100-200K standard; maintain records
- e. Repeatability - a standard in the 300K to 800K range must have a coefficient of variation (CV) of 5% or less on 10 replicates (**Refer to Operating Manual**); maintain records
- f. **THESE PARAMETERS MUST BE ACHIEVED BEFORE PROCEEDING**

### 19. Testing Samples

- a. Heat samples to 37-42°C (using a temperature control) and test within 30 min of reaching temperature
- b. Samples must not be re-used and must be discarded after use
- c. Invert 10 times, pipet 500 µL into intake chamber within 3 min
- d. Record number of cells counted for each sample

### 20. With Continuous Operation:

- a. Perform a zero check (item 17.d)
- b. Test a standard or optionally a control/sample (item 17.d) in the 600 to 800K range hourly in triplicate and determine the average, must be within 5% of the original established instrument average value (optionally, within 10% of original DMSCC average)
- c. Maintain records

### 21. Routine maintenance

- a. Maintain records

## REPORTING

### 22. Computing and Reporting Counts

- a. Count obtained x 1000 is the cell count/mL milk
- b. In reporting electronic somatic cell counts (ESCC/mL); record only first two left hand digits, raising second digit to next higher number when third digit is 6 or more
- c. Report the two left hand digits (rounded)
  - 1. If the third digit is 5 the second digit is rounded by the following rule
    - a. When the second digit is odd round up, raise the second digit by 1 (odd up, 235 to 240)
    - b. When the second digit is even round down, delete the 5 and report the second digit as is (even down, 225 to 220)
- d. If count on instrument is <100 report as <100,000 ESCC/mL
- e. If goat or camel milk is over the regulatory limit, follow confirmation procedure in PMO