

APPENDIX N BULK MILK TANKER SCREENING TEST FORM

**CHARM® ROSA SULF TEST
(Raw Commingled Cow Milk)
IMS #9-C18**

[Unless otherwise stated all tolerances are ±5%]

GENERAL REQUIREMENTS

- 1. See Appendix N General Requirements (App. N GR) items 1-8 & 15** _____

SAMPLES

- 2. See App. N GR item 9** _____

APPARATUS & REAGENTS

- 3. Equipment** _____

- a. Charm Sciences Strip Incubator:
56±1°C 8 min timer – same incubator as SL beta-lactam test;
56±1°C Charm EZ or EZ Protect+ _____

1. Clean, properly maintained and located on a level surface _____

2. Check temperature daily (day of use); maintain records _____

- a. Charm EZ and EZ Protect+ printout acceptable for daily temperature
check (annual accuracy check required); maintain records _____

3. Temperature measuring device for each incubator (App N. GR item 3) _____

4. Lid when not running tests _____

- a. 4 strip incubator lid closed (slightly sprung so that timer not active) _____

b. 2 strip incubator lid may be open _____

c. Charm EZ- reader lid closed _____

d. Charm EZ-Protect+ not applicable _____

5. Incubator Temperature: _____

6. Timer if not included in incubator
Incubation Time of internal timer: _____

- b. ROSA® Reader, ROSA Pearl Reader (with or without ROSA Barcode option),
Charm EZ, EZ-Protect+, or Charm Sciences equivalent with print out or
download of data; manual available _____

Serial Number: _____

1. ROSA Reader V1.03 or higher

a. Calibrators - 2 line for ROSA SULF Performance on SL Beta-lactam channel acceptable

Range(s):

Result

Low: _____

High: _____

b. Maintain records

2. ROSA Pearl Reader V3.00 or higher

a. Calibrators - Low and High for use in all assay channels

Range(s)

Result

Low: _____
(darker magenta)

High: _____
(lighter pink)

b. Maintain records

3. Charm EZ reader

a. Calibrators – Low and High for use in all assay channels

Range(s)

Result(s)

Low: _____
(darker magenta)

High: _____
(lighter pink)

b. Maintain records

4. Charm EZ-Protect+ each optics module in base unit

a. Calibrators – Low and High for use in all assay channels. Shall be used in all optics modules attached to EZ Protect+ base with matching serial number

b. Low and High Calibrators results – IN RANGE=1, OUT OF RANGE=0 _____
Result(s)

Low: _____
(darker magenta)

High: _____
(lighter pink)

c. Maintain records _____

5. Calibrator serial numbers match reader SN or base of EZ Protect+ _____

6. **Do not proceed if out of range.** Manufacturer should be contacted for corrective actions _____

7. Printer or computer link for hardcopy download _____

c. Pipettors - 300 µL and disposable tips (see App. N GR item 7) _____

1. **FOR SCREENING ONLY** - Single use 300 µL ROSA-pipet with overflow bulb to accurately measure amount of sample, supplied by manufacturer _____

4. Test Strips _____

a. Test Strips _____

Lot #: _____ Exp. Date: _____

QC Date: _____ By: _____

5. Sample and control agitation _____

a. Mix milk sample(s)/control(s) 25 times in 7 sec with a 1 ft movement or vortex for 10 sec at maximum setting; use within 3 min. (samples/controls must be in appropriate containers to allow the use of vortexing) _____

6. Reagent Stability and Preparation _____

a. Charm ROSA- SULF reagents must be received within 72 hours if shipped non-refrigerated; over 72 hours must be refrigerated. _____

b. Store reagents at 0.0-4.5°C, desiccant blue, maintain no longer than manufacturers expiration date. _____

1. **Do not use if desiccant indicator is white, pink or purple** _____

c. Positive Control _____

1. Manufacturer supplied, do not use after manufacturer's expiration date _____

2. Lyophilized or tablet 10 ppb Sulfadimethoxine _____

Lot # _____ Exp. Date: _____

- 3. Reconstitute 1 tablet with 5 mL Negative Control (raw milk) item 6.d _____
- 4. Tested +400 or more positive, use within 48 hours when maintained at 0.0-4.5°C _____

Lab Prep Date: _____ Lab Exp. Date: _____

- 5. Or, aliquot within 24 hours and freeze at -15°C or colder in a non-frost-free freezer or in an insulated foam container in a frost-free freezer; use within 1 month _____

Lab Prep. Date: _____ Lab Exp. Date: _____

- a. Thaw slowly overnight in refrigerator or more rapidly in cold water. Mix well until sample is homogenous _____

- 1. **Do not use if there is visible protein precipitation** _____

- b. Store at 0.0-4.5°C and use within 24 hours; do not refreeze _____

- 6. Day of use, must produce +400 or greater reading; maintain records _____

Test Value: _____

Do not proceed if out of range _____

- d. Negative Control _____

- 1. Previously tested sulfonamide negative raw milk _____
 - 2. Milk can be screened (previously tested) by the testing location making and or using the controls _____
 - 3. Negative control must test -600 or more negative _____

Sample ID: _____ Test Value: _____

Date tested: _____

- 4. Use within 72 hours when maintained at 0.0-4.5°C _____

- 5. Or, aliquot within 24 hours and freeze at -15°C or colder in a non-frost-free freezer or in an insulated foam container in a frost-free freezer; use within 2 months _____

Lab Prep. Date: _____ Lab Exp. Date: _____

- a. Thaw slowly overnight in refrigerator or more rapidly in cold water. Mix well until sample is homogeneous _____

- 1. **Do not use if there is visible protein precipitation** _____
- b. Store at 0.0-4.5°C and use within 24 hours; do not refreeze _____
- 6. Day of use must produce -600 or more negative; maintain records _____

Do not proceed if out of range _____

- e. Optional Use of Negative Control Qualifier (NCQ) – Negative Control Qualifier is used if no previously tested Negative Control (qualified raw milk, item 6.d) is available _____

- 1. Rehydrate NCQ _____
 - a. Allow bottle to warm at room temperature for 10 minutes _____
 - b. Reconstitute with 5.0 mL of 45°C potable water. Shake well _____
 - c. Allow to stand refrigerated or on ice for 15 minutes _____
 - d. Shake before use _____
 - e. Use within 48 hours when maintained at 0.0-4.5°C _____
 - f. DO NOT use NCQ to rehydrate Positive Control _____

- 2. Qualification of raw milk as Negative Control _____
 - a. Test prepared NCQ twice _____
 - 1. Both NCQ readings must be -600 or more negative _____
 - b. Test raw commingled milk twice _____
 - 1. Both raw commingled milk readings must be -600 or more negative _____
 - c. Average the two NCQ readings and the two raw commingled milk readings _____
 - 1. The average of the two commingled milk readings must be within ± 1300 of the average of the two NCQ readings _____
 - 2. If these conditions are met, the raw commingled milk is qualified as a Negative Control, item 6.d _____
 - d. Proceed to item 7 to run Performance Monitoring _____

TECHNIQUE

7. Daily Performance and Operation Check _____

- a. See App. N GR items 10.b-d. If a reader/optics module is used to test multiple tests, one of those tests shall be performance checked daily _____
- b. If using ROSA reader Versions 1.05 and higher, or ROSA-Pearl, use ESC 5 reader function to enter performance monitoring mode of reader; if using Charm EZ, use Menu to enter Performance Monitoring mode and “Perf Mon” to enter daily performance check; if using Charm EZ Protect+, use “PM” button and select “Performance Monitoring” to enter daily performance check. Each EZ Protect+ optics module in use must complete Daily Performance and Operation Check, item 7; refer to manual for directions _____
- c. Check Calibrators; items 3.b.1, 3.b.2, 3.b.3, or 3.b.4 _____
- d. Positive and negative controls must give appropriate readings prior to any sample analysis (see App. N GR item 10.b) _____
- e. Controls in-range when in performance monitoring mode, ROSA reader version 1.05 and higher, ROSA Pearl or Charm EZ or EZ Protect+ _____
- f. **Do not proceed if out of range** _____

8. Test Procedure _____

- a. Set out required number of test strips and place them in a dry labeled container at room temperature, or take out strips as needed. _____
 - 1. Discard unused test strips at the end of the day _____
- b. Label test strips, one for each test sample and each control. Avoid crushing sample compartment(s). Re-shape dented sample compartments to fit into incubator _____
- c. Mix milk sample(s)/control(s) (See Item 5 a) _____
- d. Place strip into incubator _____
 - 1. EZ reader and EZ Protect+ in incubate and read mode displays appropriate test name and message when in correct temperature range, 55-57°C _____
 - a. EZ reader displays “add milk to strip and close door”. Follow item 8.i _____
 - b. EZ Protect+ displays “add milk to strip”. Follow item 8.i within 1 minute of message display _____
- e. While holding strip flat, peel back plastic (to ‘peel to here’ line) to expose sample pad compartment. Avoid lifting the wick and sponge under tape _____
 - 1. For multiple samples, complete items 8.d-g for each sample/control, before starting test of next sample _____
 - 2. Complete all samples within 2 min of placing first strip in incubator _____

- f. Add 300 μ L of mixed sample/control to corresponding strip _____
- 1. Using pipettor (item 3.c) with new tip for each sample/control, draw up 300 μ L avoiding foam or bubbles _____
 - a. Remove tip from liquid _____
 - b. While holding the pipettor vertically, expel test portion slowly into either side well of appropriate strip _____
- 2. **FOR SCREENING ONLY** - Using new manufacturer-provided ROSA-pipet (item 3.c.1) for each sample/control _____
 - a. Squeeze top bulb while holding vertically with bulb and overflow reservoir side pointing down, draw up test portion avoiding foam and bubbles. Sample should completely fill pipet shaft and overflow into the bottom half of the overflow reservoir _____
 - b. Remove tip from liquid _____
 - c. While holding the ROSA-pipet vertically, expel test portion slowly into either side well of appropriate strip. Excess portion should remain in reservoir _____
- g. Re-seal plastic firmly around sample pad compartment. DO NOT press on compartment or sponge in compartment _____
- h. ROSA Reader and Charm EZ and EZ Protect+ (read only mode) _____
 - 1. Close lid and latch ROSA incubator to start automatic timer in the incubator. If no automatic timer in incubator, set external timer for 8 min. Incubate 8 min not to exceed 9 min. _____
 - 2. At end of incubation visually inspect C (Control) line an absent C line, a partial C line or an indistinct C line indicates an invalid test; and the sample/control must be re-tested _____
 - 3. Insert only valid test(s) in reader _____
 - a. ROSA reader set to appropriate channel _____
 - 1. SULF slow blink for ROSA SULF test _____
 - 2. Press ENTER, reading and interpretation appear in 5 sec, read strips within 5 min of completion of incubation _____
 - b. Charm EZ and EZ Protect+ automatically sets channel when color/bar coded strip inserted _____
 - 1. Charm EZ- Close door; reading and interpretation appear in 5 sec, read strips within 5 min of completion of incubation _____

2. Charm EZ Protect+ - reading and interpretation appear in 5 sec on strip insertion, detection and entering required information _____

i. Charm EZ and EZ Protect+ (incubate and read mode) _____

1. Charm EZ and EZ Protect+ automatically sets channel and incubator temperature when color/bar coded strip inserted. Optionally enter sample ID and/or operator. Wait for "Add milk to strip" message indicating incubator temperature is in range _____
2. Peel strip (8.e), add milk (8.f) and re-seal strip (8.g) _____
3. Start Reading _____
 - a. Charm EZ- Close door to begin _____
 - b. EZ-Protect+- Count starts automatically on strip and milk flow detection and times out if milk not detected in specified time _____
4. Charm EZ and EZ Protect+ automatically prompts for further testing when positive _____

9. Interpretation with Reader _____

- a. If there is a negative or zero reading on the reader, sample is a Negative (NF) _____
- b. If there is a positive reading on the reader, sample is an Initial Positive _____

10. Verification of Initial Positive Tanker Samples Done at the Same Testing Facility (see App. N GR item 11) _____

- a. Set out four test strips and label as negative control, positive control, and two strips with the initial positive sample ID. Avoid crushing sample compartment(s) _____
- b. Mix milk sample(s)/control(s) (See Item 5.a) _____
- c. Place a labeled strip into incubator _____
- d. While holding strip flat, peel back plastic (to 'peel to here' line) to expose sample pad compartment. Avoid lifting the wick and sponge under tape _____
 1. Complete items 10.c-f for each sample/control, before starting test of next sample/control _____
 2. Complete pipetting all samples within 2 min of placing first strip in incubator _____
- e. Add 300 µL of mixed sample/control to corresponding strip _____

1. Using pipettor (item 3.c) with new tip for each sample/control, draw up 300 µL avoiding foam or bubbles _____
 - a. Remove tip from liquid _____
 - b. While holding the pipettor vertically, expel test portion slowly into either side well of appropriate strip _____
- f. Re-seal plastic firmly around sample pad compartment _____
- g. ROSA Reader, Charm EZ and EZ Protect+ (must confirm in read only mode) _____
 1. Close lid and latch ROSA incubator to start automatic timer in the incubator. If no automatic timer in incubator, set external timer for 8 min. Incubate 8 min not to exceed 9 min _____
 2. At end of incubation remove strips from incubator. Visually inspect each strip C (Control) line. An absent C line, a partial C line or an indistinct C line indicates an invalid test; and the sample/control must be re-tested _____
 3. Insert only valid test(s) in reader _____
 - a. ROSA reader set to appropriate channel _____
 1. SULF slow blink for ROSA SULF test _____
 2. Press ENTER, reading and interpretation appear in 5 sec, read strips within 5 min of completion of incubation _____
 3. Positive and negative controls must give appropriate readings prior to any sample analysis (see items 6. c-d) _____
 - a. Controls in-range when in performance monitoring mode, ROSA reader version 1.05 and higher, ROSA Pearl _____
 - b. Charm EZ and EZ Protect+ automatically sets channel when color/bar coded strip inserted _____
 1. From initial positive screen, on Charm EZ press Continue or Confirm Later on Charm EZ Protect+ press Verify or Verify Later _____
 - a. If press Confirm/Verify Later, then call the confirmation routine back by pressing Main Menu _____
 - b. On Charm EZ- Select Perf Mon. then from next menu select CNFM-Later button. On Charm EZ Protect+ select PM icon from main page. Select Verify/Confirm Positive _____
 - c. Select initial positive sample ID and click ok _____

2. Follow Prompts to insert negative and then positive control. Click Continue Confirmation upon successful completion of positive control. _____
3. Follow prompts to enter the first duplicate sample. Insert the strip and on ROSA and Charm EZ readers then close door; reading and interpretation appear in 5 sec. Follow prompts to enter the second reading and interpretation appear in 5 sec. Read strips within 5 min of completion of incubation _____
4. Charm EZ and Charm EZ Protect+ automatically determines Presumptive Positive or Not Found. _____

11. **Confirmation of Presumptive Positive Tanker Samples (see App. N GR item 12) [Only in an accredited laboratory or by a CIS]** _____
12. **Traceback of Producer(s) on a Confirmed Positive Tanker (see App. N GR item 13) [Only in an accredited laboratory or by a CIS (refer to M-a-85 current revision for a listing of test kits to assure equivalence)]** _____
13. **Re-instatement of Producer(s) [Only in an accredited laboratory or by a CIS (refer to M-a-85 current revision for a listing of test kits to assure equivalence)]** _____
14. **Reporting (see App. N GR item 14)** _____